

Bone Cement

Knochenzement

190748

What are the contents of the package?

The package contents of Solid (50 g - 45 g Polymethylmethacrylate + 4,6 g Zirconium dioxide + 0,4 g Activable peroxyde benzoyl) and Liquid (ampoul of 21 ml - 17,7 ml Methyl methacrylate + 2,9 ml n-butyl methacrylate + 0,4 ml N,N-dimethylparatoluidine).

How is it packed and how often can you use it?

Sterile packed and for single use.

Was ist in der Packung enthalten?

Es sind Feststoffe (50 g - 45 g Polymethylmethacrylat + 4,6 g Zirkoniumdioxid + 0,4 g aktivierbare Benzoyl-peroxyde) und Flüssigstoffe (Ampulle von 21 ml - 17,7 ml Methylmethacrylat + 2,9 ml n-Butyl-Methacrylat + 0,4 ml N,N-Dimethylparatoluidin) enthalten.

Wie ist er verpackt und wie oft kann man ihn verwenden?

Er ist steril verpackt und nur einmal verwendbar.

EICKEMEYER®

GERMANY | EltstraÙe 8 | 78532 Tuttlingen | Deutschland | T +49 7461 96 580 0 | F +49 7461 96 580 90 | info@eickemeyer.de | www.eickemeyer.de

SWITZERLAND | Sandgrube 29 | 9050 Appenzell | Schweiz | T 071 788 23 13 | F 071 788 23 14 | info@eickemeyer.ch | www.eickemeyer.ch

UNITED KINGDOM | 3 Windmill Business Village | Brooklands Close | Sunbury-on-Thames | Surrey | TW16 7DY | UK | T 020 8891 2007 | www.eickemeyer.co.uk

DENMARK | Solbakken 26 | Hammelev | 6500 Vojens | Danmark | T +45 70 20 50 19 | info@eickemeyer.dk | www.eickemeyer.dk

NETHERLANDS | Bellweg 44 | 4104 BJ Culemborg | Nederland | T +31 345 58 9400 | info@eickemeyer.nl | www.eickemeyer.nl

ITALY | Via G. Verdi 8 | 65015 Montesilvano (PE) | Italia | T 085 935 4078 | F 085 935 9471 | info@eickemeyer.it | www.eickemeyer.it

POLAND | Al. Jana Pawła II 27 | 00-867 Warszawa | Polska | T 22 185 55 76 | F 22 185 59 40 | info@eickemeyer.pl | www.eickemeyer.pl

CANADA | 617 Douro Street, Suite #205 | Stratford | Ont. | Canada | N5A 0B5 | T +1 519 273 5558 | F +1 519 271 7114 | info@eickemeyervet.ca | www.eickemeyercanada.ca